

## Layout



## Switch Control

It is possible to connect a switch to the Omni  
The switch should be connected to the port on the bottom left hand corner.



It is important that the switch is placed in the correct position to give easy access and control of the system.

**Care must be taken when connecting the switch, do not force the switch connection as this will cause damage to the control unit.**

## Charging Omni

**Please ONLY use the supplied charger** – please take care when connecting & disconnecting the charge lead and don't force the connection, **using force will cause damage to the charge port**

When charging, the LED of battery status will remain red. When the battery is fully charged, the LED will be green. When the LED is green, you can disconnect the battery charger.

The full charging time of the battery can vary, in most situations the unit will require charging daily.

**If the device has been allowed to go completely flat, you might need to charge the device for 5-10 minutes before it can be used again. You will also need to turn the unit back on using the stylus whilst the charge is connected for the unit to start charging correctly.**

## Turning Omni off & on

### Turning off

Using the provided stylus press **gently** into the 2<sup>nd</sup> hole down from the top on the right-hand side of the unit for 3 seconds.



### Turning On

Using the provided stylus **press gently** into the 2<sup>nd</sup> hole down from the top on the right-hand side of the unit for 3 seconds.

After a few moments you will get a prompt on screen of 'Start menu' or 'Setup' please select 'Start Menu'.

## Controlling Infrared

Omni can control devices like Televisions, DVD players & other devices that operate via an infrared remote control. These devices need to be in the same room as the Omni unit.

If you are having difficulty in controlling these devices please check:

- the Omni is in the same room
- and that it is pointing roughly towards the device you are trying to control
- the Omni is charged.



Tel: 0113 2564572  
Email: [info@dh2solutions.co.uk](mailto:info@dh2solutions.co.uk)  
Web: [www.dh2solutions.co.uk](http://www.dh2solutions.co.uk)

## Omni User Guide

# Rebooting the omni

If you are experiencing any difficulties with the omni the first step would be to turn the Omni off and back on again.