

Name: Floorstand PentaLock Tele (no weight)

Model Number: #17.1100 / #17.1101

Instructions ID: #00.5087~C


## Floorstand PentaLock Tele



Rehadapt



# ASSEMBLY INSTRUCTIONS




Replacement / Austausch




## Mounting kit / Montageset

- |                                 |                                        |
|---------------------------------|----------------------------------------|
| (A) = 1x FS base<br>FS Basis    | (E) = 1x FS arm<br>FS Arm              |
| (B) = 5x FS caster<br>FS Rollen | (F) = 1x T-handle<br>T-Griff           |
| (C) = 1x UDS<br>UGA             | (G) = 1x M8x50 screw<br>M8x50 Schraube |
| (D) = 1x FS column<br>FS Säule  | (H) = 1x M8x35 screw<br>M8x35 Schraube |


1



2









3



## Attention / Achtung



# ASSEMBLY INSTRUCTIONS





# OPERATING INSTRUCTIONS



## Attention / Achtung



Do not adjust/operate any hinge with a device mounted!

Remove the device before adjusting!



Betätigung/Verstellung jeglicher Gelenke nur mit abgenommenem Gerät!

Bitte Gerät vorher abnehmen!




## Information / Hinweis



### Telescopic column

The height of the column is adjustable so the device can be placed precisely at the desired position.



### Teleskopierbare Säule

In der Höhe teleskopierbare Säule zum Positionieren des Gerätes in gewünschte Position.